

## Shade Indicator Prismatik ThinPress Ingots

ThinPress Shades	TP1	TP2	TP3	TP4	TP5
Thickness	0.3 - 0.5	Thick or Thin	0.3 - 0.5	Thick or Thin	0.3 - 0.5
Corresponding Shades	<b>Bleach</b> ex: 010, OM1, OM2, OM3, NW0, NW0.5, B0, i1-i8	<b>Bleach</b> ex: A0, 020, 030, NW0, NW0.5, B0, OM2, OM3, XL	A1, A2	A3, A3.5, A4	B1, B2

ThinPress Shades	TP6	TP7	TP8	TP9	TP10
Thickness	Thick or Thin	Thick or Thin	Thick or Thin	Thick or Thin	Thick or Thin
Corresponding Shades	C1, C2, 040	D2, D3	A1, A2, B1, B2	B3, B4	C3, C4, D4

\* Thin used for pressings as thin as 0.3mm - 0.5mm, no prep or minimal prep.  
(Note - TP1, TP3, TP5, are used to maximize opacity and block out dark underlying tooth shades)

\* Thick used for pressings thicker than 0.6mm, ideal veneer prep or full crown prep.  
(Note - TP2, TP4, TP6, TP7, TP8, TP9, TP10 contain more translucency)

## Correction Powders

Correction Body Shades	CP1	CP2	CP3	CP4	CP5	CP6	CP7
Corresponding Shades	Bleach	A1, B1, A2, B2	A3, A3.5, A4	C1, C2	D2, D3	B3, B4	C3, C4, D4