

Stop snoring. Sleep better. Get aveoTSD®.

aveoTSD®

anti-snoring aid

Good Health Through Quality Sleep

**The brilliantly
simple solution for
problem snoring.**

aveoTSD®




Patients:

Please print this sheet and take it to your dentist or physician who can prescribe you an aveoTSD.


Dentists/Physicians:

*Call 800-334-1979 to order.
Visit getaveo.com for more info.*

How the aveoTSD works



In this image, the tongue falls into the back of the airway as a person sleeps. This blocks the airway, leading to snoring.



This image shows the aveoTSD holding the tongue gently forward, preventing it from falling back and obstructing the airway. Note how the airway is now open and clear. This stops or greatly reduces snoring.



**GLIDEWELL
LABORATORIES**

"The Standard of Care"

800-334-1979 • getaveo.com

GL-2267-0710