



APPOINTMENT #2: MULTI-UNIT ABUTMENTS, IMPLANT VERIFICATION JIG, ABUTMENT-LEVEL IMPRESSIONS, WAX RIM AND BITE

What you will receive:

- Master cast with multi-unit abutments
- Implant verification jig (IVJ)
- Custom impression tray
- Wax rim with six anterior teeth
- Prosthetic screws
- Opposing cast

What you will need:

PART A: Multi-Unit Abutments, IVJ and Final Impression

- Implant prosthetic driver
- Multi-unit abutment driver
- Torque wrench
- X-ray bitewing tab
- Flowable composite
- Curing light
- Carbide bur
- Light- and heavy-body VPS
- VPS adhesive
- Periphery wax

PART B: Wax Rim

- Implant prosthetic driver
- Burner (heat source)
- Buffalo knife
- Occlusal plane former
- Fox plane
- Wooden tongue depressor
- No. 7 wax spatula
- Baseplate wax
- Lab handle and scalpel
- Bite registration (rigid impression material)
- Camera

What you will send back to Glidewell:

- Final impression with the IVJ
- Adjusted wax rim
- Bite registration
- Opposing cast
- Smile photo of the patient wearing the adjusted wax rim
- List of additional changes for setup

What you will receive next:

- Mounted master cast and opposing cast
- Wax setup
- Prosthetic screws